

PRODUCT NON-CONFORMITY RECORD (CAPA)		YANTRAM MEDTECH PVT LTD
Doc. No. YMPL-PNC-CAPA-F01 REV NO: 00		Page 1 of 1

Department	Quality Control	Date	
Details of Non- Conformity			
NV_CORE: Charging LED is not turn On while device charging.			
Observed By	M V R Amar		
Action taken to correct (remedial) non-conformity			
Replaced the Charging LED.			
Root Cause Analysis (main reason)			
Led damaged.			
Corrective Action Taken to avoid the recurrence			
Check the any solder bridges / solder balls / device drops from height.			
Responsibility	M V R Amar	Target Date	
Preventive Action Taken to avoid the occurrence			
Maintain proper soldering/ handle the device with care.			
Remarks regarding effectiveness & Sign			
Replaced the Charging LED and Tested.			